

Considerations when analysing PFAS containing samples

Matt James, Tony Edge

Agenda

01

- Introduction to PFAS

02

- Contamination issues with PFAS

03

- Use of delay columns

04

- Sample preparation techniques for environmental analysis

05

- Conclusions

06

Introduction to PFAS

- Per- and polyfluoroalkyl substances are a group of synthetic organofluorine chemical compounds that have multiple fluorine atoms attached to an alkyl chain
- Often referred to as forever chemicals due to their longevity in the environment, they have also been found in a wide range of marine life
- US EPA toxicity database, DSSTox, lists 14,735 unique PFAS chemical compounds,^[7] while PubChem lists approximately 6 million
- Many PFASs were used for their enhanced water-resistant properties
- Residues have been detected in humans and wildlife, prompting concern about impacts to health.^{[9][10][11]} According to the National Academies of Sciences, Engineering, and Medicine
- PFAS exposure is linked to increased risk of dyslipidemia (abnormally high cholesterol), suboptimal antibody response, reduced infant and fetal growth, and higher rates of kidney cancer.^[12]

Challenges to developing PFAS workflow

- PFAS are ubiquitous within the laboratory environment
 - Very challenging to prevent environmental contamination
 - Use of delay columns reduces potential interferences derived from analytical system
- Compounds are thought to be highly toxic and so detection limits are very low
 - Typically 1 ng/L
 - Often requires pre-concentration using SPE, e.g. factor of 250x according to EPA 537.1
 - Low detection limits require LC-MS/MS
- Some PFAS compounds will adsorb to plastic and glass surfaces, reducing recovery
 - PFAS compounds can also stick to Aluminium
- Toxicity concerns mean that extra handling precautions should be taken
 - This can be a source of possible contamination


List of EPA 537.1 compounds and MRLs

TABLE 5. DLs AND LCMRLs IN REAGENT WATER

Analyte	Fortified Conc. (ng/L) ^a	DL ^b (ng/L)	LCMRL ^c (ng/L)
PFBS	4.0	1.8	6.3
PFHxA	4.0	1.0	1.7
HFPO-DA	4.0	1.9	4.3
PFHpA	4.0	0.71	0.63
PFHxS	4.0	1.4	2.4
ADONA	4.0	0.88	0.55
PFOA	4.0	0.53	0.82
PFOS	4.0	1.1	2.7
PFNA	4.0	0.70	0.83
9Cl-PF3ONS	4.0	1.4	1.8
PFDA	4.0	1.6	3.3
NMeFOSAA	4.0	2.4	4.3
PFUnA	4.0	1.6	5.2
NEtFOSAA	4.0	2.8	4.8
11Cl-PF3OUdS	4.0	1.5	1.5
PFDoA	4.0	1.2	1.3
PFTTrDA	4.0	0.72	0.53
PFTA	4.0	1.1	1.2


PFAS Method: contamination testing

- EPA 537.1 method
- This method is used for contamination testing



Column:	Avantor® ACE® UltraCore SuperC18, 2.5 µm, 100 x 2.1 mm														
Delay column:	Avantor® ACE® PFAS Delay Column, 50 x 2.1 mm														
Mobile Phases:	A: 20 mM ammonium acetate B: MeOH														
	<table><tr><th>Time (mins)</th><th>% B</th></tr><tr><td>0</td><td>5</td></tr><tr><td>0.1</td><td>40</td></tr><tr><td>8.5</td><td>95</td></tr><tr><td>10.5</td><td>95</td></tr><tr><td>10.6</td><td>5</td></tr><tr><td>13.2</td><td>5</td></tr></table>	Time (mins)	% B	0	5	0.1	40	8.5	95	10.5	95	10.6	5	13.2	5
Time (mins)	% B														
0	5														
0.1	40														
8.5	95														
10.5	95														
10.6	5														
13.2	5														
Flow Rate:	0.4 mL/min														
Temperature:	40 °C														
Detection:	Sciex QTRAP® 6500+ LC-MS/MS Ionisation mode: ESI, negative mode; Source temperature: 450 °C; Curtain gas: 30 psig; Ionspray™ source voltage: -4500 V														
Sample:	EPA 537.1 mix (18 analytes), 3 internal standards and 4 surrogates														


PFAS Method: Delay column use



Contaminated LC-MS system

Blank run

No delay column



Contaminated LC-MS system


Injection of 0.5 ng/mL PFAS standards

No delay column

- Sample PFAS coelute with system background PFAS that accumulate on-column during post gradient re-equilibration
- Results in inaccurate quantitation

PFAS Method: Delay column use - PFOA

XIC from Data_5.mill (sample 1) - Delay column A/N data_Analytical only, no injection_005, MRM(18 transitions): PFOA_1 (412.8/368.9), Gaussian smoothed (1.0 points)




Contaminated LC-MS system

Blank run

No delay column

XIC from Data_9.mill (sample 1) - Delay column A/N data_Analytical only, PFAS 0.5 ng/mL_009, MRM(18 transitions): PFOA_1 (412.8/368.9), Gaussian smoothed (1.0 points)




Contaminated LC-MS system

Injection of 0.5 ng/mL PFAS standards

No delay column

XIC from Data_14.mill (sample 1) - Delay column A/N data_Delay 2209/39, no injection_014, MRM(18 transitions): PFOA_1 (412.8/368.9), Gaussian smoothed (1.0 points)




Contaminated LC-MS system

Blank run

Delay column installed

XIC from Data_17.mill (sample 1) - Delay column A/N data_Delay 2209/39, PFAS 0.5 ng/mL_017, MRM(18 transitions): PFOA_1 (412.8/368.9), Gaussian smoothed (1.0 points)




Contaminated LC-MS system

Injection of 0.5 ng/mL PFAS standards

Delay column installed


PFAS Method: Delay column use



Contaminated LC-MS system

Blank run

Delay column installed



Contaminated LC-MS system



Injection of 0.5 ng/mL PFAS standards

Delay column installed

- Use of a delay column retains system background PFAS until after sample PFAS have eluted


Tracking system contamination: solvent tubing

3.0 cm of solvent tubing extracted into 1.5 mL MeOH




- Similar PFAS profile observed to that obtained when no delay column was used.
- Solvent tubing is a significant source of system related PFAS.
- Switch tubing to PEEK

Tracking system contamination: Blanks




10 μ L injection of 1 mL MeOH in transfer tube



Blank – no injection


- Blanks are relatively clean.
- Conclusion – solvent tubing contributes significantly to system background PFAS levels

Sample Prep: N₂ Blowdown



- 8 mL MeOH blown down to dryness in 10 mL polypropylene transfer tube.
- N₂, 40L/min @ 60 °C
- Reconstituted in 1 mL MeOH
- All samples contaminate with PFHxA
- Level of contamination is consistent between positions
- Level of contamination consistent over time
- Originates from Blowdown rig or N₂ source?


Sample Prep: N₂ Blowdown



- Blowdown rig was removed.
- N₂ direct from N₂ generator via polypropylene tubing and glass pipette
- No heat, so longer required to blow down, therefore increased signal!
- N₂ source appears to be the culprit.
- N₂ generator contains PTFE filters – potentially the source?
- Use activated carbon inline filter?

Sample Prep: Vials & caps matter!

- ADONA was observed in blank MeOH injections



Use PP vials & caps!

Polyimide cap is source

Cap or vial?


ALS OK

system OK

Testing sequence


Other sources of contamination

- Samples of common laboratory consumables were extracted in 1 mL MeOH and analysed




Syringe Filters

- Many samples contain particulates
 - Sub-micron particles → → → → large particles/micro-organisms
- Can damage LC system components & analytical column
 - Blockages
 - Peak distortion
 - Increase in back pressure
 - Reduced column lifetime
- Syringe filters provide fast, cheap & convenient solution




Filter quality

- Not all syringe filters are created equal!
- Important to use filters with high extraction efficiency





Extraction Efficiency Test: J.T.Baker 0.22 µm Filters



Extraction efficiencies determined by filtering a 0.01% latex solutions polystyrene beads through syringe filters by spectrophotometric analysis (UV, 272 nm). Triplicate analyses performed for multiple batches.

Filter quality: Extractables

- Poor Quality filters can contaminate samples
 - Ghost peaks
 - Interfere with target analyte quantification
- Using high quality syringe filters protects sample integrity.
- Filters easily tested for extractables by LC-UV:




Method:

1 mL of H₂O or MeOH passed through filter.
Eluent analysed by HPLC

Column: Avantor® ACE® 3 C18, 150 x 4.6 mm
Mobile phase A: H₂O
Mobile phase B: MeCN
Gradient: 5 to 100% B in 15 mins.
Flow rate: 1 mL/min
Temperature: 30 °C
Detection: UV, 214 nm
Injection volume: 100 µL

SPE - Visualising the process (schematically):



Solid Phase Extraction (SPE)

- Method development required
 - Relatively easy to do
 - Many manual manipulation steps
- Sample preconcentration
 - Possible to preconcentrate >100 times
 - Method development requires some skill
 - Solvent usage less than liquid-liquid extraction
- Minimal issues with matrix as removes more of matrix than other techniques
- Applicable to a wide range of compounds with differing log D
- Excellent recoveries



Selecting a sorbent

Interaction Mode	Interaction Energy (kcal/mol)
Ionic Interactions	50 – 200
Polar Interactions	3 – 10
Non-polar Interactions	1 – 5

- Stronger retention offered by ion exchange
- Combination of interactions provides best for selectivity

SPE – Considerations for the Five-Stage Process

1. Condition

- Strongly eluting solvent
- Suitable volume

2. Equilibrate

- Weakly eluting solvent
- Suitable volume

3. Load

- Ionization state
- Concentration of analyte vs bed weight
- Cleanliness of matrix

4. Wash(es)

- Elutropic strength
- Suitable volume

5. Elute

- Elutropic strength
- Suitable volume

Conclusion

01

- Introduction to PFAS

02

- Contamination issues with PFAS

03

- Use of delay columns

04

- Sample preparation techniques for environmental analysis

05

- Conclusions

06

- Acknowledgements

Acknowledgements

- Avantor colleagues: Gemma Lo, Arianne Soliven, Gemma Howse and Amelia Knight

Thank you

Anthony.edge@avantorscience.com

chromsupport@avantorsciences.com